

Microprocessor LABORATORY CENTRIFUGE MACHINE



High Speed • Digital Display
• Heavy Duty Body

Specifications

- ✓ **RPM Range:** 3000 – 6000 RPM
- ✓ **Capacity:** 8 x 15 ml Tubes
16 x 15 ml Tubes
24 x 15 ml Tubes
- ✓ **Display:** Microprocessor Based Digital Display
- ✓ **Timer:** 0 – 60 Minutes Adjustable
- ✓ **Motor:** Heavy Duty Copper Winding Motor
- ✓ **Body:** Powder Coated Heavy Metal Body
- ✓ **Rotor:** Aluminium Rotor with Tube Holders
- ✓ **Operation:** Low Noise & Low Vibration
- ✓ **Safety:** Lid Safety Lock System
- ✓ **Control:** Start / Stop Button with RPM Setting
- ✓ **Power Supply:** 220V AC, 50 Hz
- ✓ **Application:** Blood / Serum Separation for Pathology Labs & Clinics
- ✓ **Compact Table Top Design**
- ✓ **Energy Efficient**
- ✓ **Long Life Motor**
- ✓ **Easy Maintenance**



Marketed By:
