

UV-VIZ Basic TRANSILLUMINATOR



Make : ORANGE

Model : UV-VIZ Basic Transilluminator

Cat.No : UVG3657

UV Source : 365nm Wavelength & Visible light

Viewing area : 20x14cm / 20x14cm

UV Tube : 4x8 watts UV Tube & Visible Light

Housing : Anti Corrosion treated with Powder coated Mild Steel

UV Safety Shield : Clear Transparent lid for excellent clarity

Input : 220V AC, 50-60HZ

Fuse : 2.5Amp

Anti Slip rubber foot provided for base

Instruction Manual & Warranty card



UVG3657