



DispensMate-Pro

Premium Bottle-Top Dispenser

Designed Safe, Sturdy & Stable dispensing directly from the bottle in the range of volumes up to 100 mL

- ① Fast, precise, and reproducible volume locking mechanism on a clear scale
- ② Reagent recovery function reduces reagent waste and prevents dripping when not in use
- ③ 6 adapters to comfortably fit most laboratory reagent bottles
- ④ Adjustable aspirating tube for different bottle heights

Features

- High chemical compatibility with use of inert materials and PTFE piston
- Easy to disassemble and clean for periodic maintenance
- Fully autoclavable at 121°C to ensure sterility
- Available in 5 dispensing volume ranges from 0.5mL to 100mL
- Supplied with most popular adapters (S40, GL32, GL38, GL25, GL28)

Specifications

Cat. No.	Volume Range	Graduation	Accuracy error		Precision error	
	mL	mL	μL	%	μL	%
7032212001	0.5–5	0.1	± 25	± 0.5	± 10	± 0.2
7032212002	1.0–10	0.2	± 50	± 0.5	± 20	± 0.2
7032212003	2.5–25	0.5	± 125	± 0.5	± 50	± 0.2
7032111004	5.0–50	1.0	± 250	± 0.5	± 100	± 0.2
7032111005	10–100	2.0	± 500	± 0.5	± 200	± 0.2

Cat. No.	17400131	17400132	17400133	17400134	17000069	17400037	17400021	17100073	17100553
Descriptions	GL32,Adapter, 45/32mm	GL38,Adapter, 45/38mm	GL25,Adapter, 32/25mm	GL28,Adapter, 32/28mm	S40, Adapter, 45/40mm	Reagent Bottle (Brown, 1L)	Filling tube (m)	Discharge tube (m)	GL45,adapter, 45mm