


GRUN MIDI Plus Horizontal Immersed Gel Electrophoresis Unit



Specification:

Make: ORANGE

Cat:GRUN50H

Dimensions: 310x150x90mm

Buffer Volume: 650ml

Number of samples: 50

Small gel tray:60x60mm (2 Nos)

Medium gel tray:120x60mm (1 No)

Large gel tray:120x120mm (1 No)

25/11 wells 1mm thick comb (2 Nos)

18/8 wells 1.5mm thick comb (1 No)

13/6 wells 1.5mm thick comb (1 No)

3/2 wells 2mm thick comb (1 No)

Buffer chamber with safety lid

Positive and negative platinum electrodes

Gel casting stand(138x138x36mm): 1 No

Connecting cord: Red and Black

Gel Scoop : 1 No

Instruction Manual & Warranty Card



GRUN50H