

Electrophoresis

Power Supply 300V/200mA



Cat.No : DPS2003

Output range : 10–300 V, fully adjustable in 1 V steps, 1–200 mA, fully adjustable in 1 mA steps

Type of output : Constant voltage or constant current

Output terminals : 4 pair recessed banana jacks in parallel

Timer : 1 min–999 min , fully adjustable

Start /Stop Function : Yes

Display : 3-digit LED

Operating conditions : 0–40°C, 0–90% humidity in the absence of condensation
overload/short-circuit detection,

Input protection : 1Amp. Fuse on both hot and neutral

Input power (actual) : 220 VAC/50hz

Instruction Manual & Warranty card



DPS2003