

PORTABLE AUTOCLAVE ALUMINIUM

TAI 901

- All Body Made of Aluminium.
- Fitted with pressure gauge.
- The unit has spring loaded safety valve for pressure adjustment and a dead weight safety valve.
- It is hydraulically tested up to 30 PSI & is supplied with Electric fitting, Cord & Plug, Supplied with tripod Stand.
- Works on 220/330 volt AC Supply.

MODEL	SIZE (MM)	SIZE (INCHES)	CAPACITY (LITER)
901	300x300	12x12	20 Liters
901	300x375	12x15	25 Liters



TAI 901 - A PORTABLE 'P' TYPE/COOKER TYPE AUTOCLAVE

- All Body Made of Aluminium, powder coated or stainless steel.
- Fitted with pressure gauge.
- The unit has spring loaded safety valve for pressure adjustment and a dead weight safety valve.
- It is hydraulically tested up to 30 PSI & is supplied with Electric fitting, Cord & Plug, Supplied with tripod Stand.
- Works on 220/330 volt AC Supply.

MODEL	SIZE	SIZE (INCHES)	CAPACITY (LITER)
901 (A)	225X275	09X11	13 Liters
901 (A)	300X300	12X12	20 Liters
901 (A)	300X500	12X20	40 Liters



ALUMINUM



POWDER COATED



STAINLESS STEEL