

GRUN Mini Plus Horizontal gel Electrophoresis Unit



Specifications:

Model # : GRUN32H

Construction material: Poly Carbonate

Gel Size (W x L) : 10 x 10cm

Buffer Required : 400ml

Platinum Electrode : 1 set anode and cathode
removable electrodes

Standard Combs : 8 & 16 well height adjustable combs
(2 Nos Each), 1mm , 1.5mm Thickness

Maximum Sample capacity: 32 Samples

Comb Slots : 5 Slots

Gel Caster : Sillicon Rubber Casting Dam

Connecting cord : 1 Pair (Red and Black) Included as
standard

